


Carbon Black Reagents for ASTM D1510

- Sodium thiosulfate, 0.0394 N
- Iodine, 0.0473 N
- Prepared from ACS grade materials
- Matrix prepared using double-deionized, ASTM Type I Water



spex.com

Phone: +1.732.549.7144 • +1.800.LAB.SPEX
Fax: +1.732.603.9647
spexsales@antylia.com

Connect with us



Spex CertiPrep is an
Antylia Scientific company.
Find out more at antylia.com.



Carbon Black Reagents for ASTM D1510

Details matter...

Our sodium thiosulfate solutions are prepared from ACS grade, micro-crystalline materials. In order to maximize shelf life, our matrix is prepared using double-deionized, ASTM Type I Water.

Our iodine solutions are prepared from ACS grade potassium iodide and crystalline elemental iodine. To guarantee a clean and stable product, our matrix is prepared using double-deionized, ASTM Type I Water.

All solutions are prepared gravimetrically using high accuracy analytical balances to ensure precise target concentrations. Each batch is thoroughly homogenized using a high speed industrial mixer to ensure reliable results from the first bottle to the last.

We are titrating our samples on our automated titrator. The automated dosing drive uses 10,000 steps over a 20 mL volume, so its dosing increment *can be* as small as 2 μL . For these applications, we are using a minimum dose of 10 μL for the sodium thiosulfate endpoint and 4 μL for the iodine endpoint. These settings achieve the extremely precise measurements for each titration while also staying within the parameters of the dosing unit.



Inorganic Certified
Reference Materials



Gravimetrically
prepared



Supplied with a
Certificate of Analysis



ISO Accredited
Standards

Carbon Black Reagents for ASTM D1510

As stated on our Certificate of Analysis, the sodium thiosulfate is run against a 0.1 N potassium dichromate solution. The exact normality of this solution is calculated by comparing it to NIST potassium dichromate. A set of 6 samples are run that must all be within the nominal value of $0.0394\text{ N} \pm 0.00008\text{ N}$. The certified sodium thiosulfate is then used to titrate iodine. A set of 3 samples are run that must all be within the nominal value of $0.0473\text{ N} \pm 0.00003\text{ N}$. Before releasing either of these reagents for packaging, we run QC checks with a previous lot to ensure accuracy over time.

Want a free sample?

Contact us today at spexsales@antylia.com or +1.732.549.7144 to request a free sample for method validation.

Carbon Black Reagents			
Description	Packaging	Volume	Part #
0.0394 N Sodium Thiosulfate	Cubitainer	1 gallon	182002
0.0473 N Iodine	Amber Glass Bottle	1 gallon	183134

Larger volume containers available upon request.

spex.com

Phone: +1.732.549.7144 • +1.800.LAB.SPEX
Fax: +1.732.603.9647
spexsales@antylia.com

4772J

Connect with us



Spex CertiPrep is an
Antylia Scientific company.
Find out more at antylia.com.

