

SPEX

Explore Reliable PFAS Standards

Working with reference
materials you can trust
is critical to success

PFAS testing can be highly challenging which is why we make it easier with our convenient portfolio of PFAS reference materials. Advance environmental and health safety with standards you can trust.

Featured Products:

- PFAS-PFOS
- PFAS-PFOA
- PFAS-PFHxS
- PFAS-PFNA

PFAS CRMs:

- 29 individual products
- 100 ppm in 95:5 Methanol:Water, Trace Sodium Hydroxide
- Tested for EPA methods 533 and 537.1
- Useful for various matrices
- Custom standards of other PFAS compounds available on demand



[SPEX.COM](https://www.speex.com)

+1.732.549.7144

+1.800.LAB.SPEX

PFAS Singles

Concentration: 100 ppm in 95:5 Methanol:Water, Trace Sodium Hydroxide

Part Number	Compounds	CAS No.	Method 533	Method 537.1
PFAS-PFEESA	Perfluoro(2-ethoxyethane)sulfonic acid	113507-87-2	x	
PFAS-HFPO-DA	hexafluoropropylene oxide dimer acid	13252-13-6	x	x
PFAS-NFDHA	Nonafluoro-3,6-dioxaheptanoic acid	151772-58-6	x	
PFAS-PFOS	Perfluorooctanesulfonic acid	1763-23-1	x	x
PFAS-PFUnA	Perfluoroundecanoic acid	2058-94-8	x	x
PFAS-NMeFOSAA	N-methyl perfluorooctanesulfonamidoacetic acid	2355-31-9		x
PFAS-PFPeA	Perfluoropentanoic acid	2706-90-3	x	
PFAS-PFPeS	Perfluoropentanesulfonic acid	2706-91-4	x	
PFAS-6:2FTS	1H,1H, 2H, 2H-Perfluorooctane sulfonic acid	27619-97-2	x	
PFAS-NEtFOSAA	N-ethyl perfluorooctanesulfonamidoacetic acid	2991-50-6	x	
PFAS-PFHxA	Perfluorohexanoic acid	307-24-4	x	x
PFAS-PFDoA	Perfluorododecanoic acid	307-55-1	x	x
PFAS-PFOA	Perfluorooctanoic acid	335-67-1	x	x
PFAS-PFDA	Perfluorodecanoic acid	335-76-2	x	x
PFAS-PFHxS	Perfluorohexanesulfonic acid	355-46-4	x	x
PFAS-PFBA	Perfluorobutanoic acid	375-22-4	x	
PFAS-PFBS	Perfluorobutanesulfonic acid	375-73-5	x	x
PFAS-PFHpA	Perfluoroheptanoic acid	375-85-9	x	x
PFAS-PFHpS	Perfluoroheptanesulfonic acid	375-92-8	x	
PFAS-PFNA	Perfluorononanoic acid	375-95-1	x	x
PFAS-PFTA	Perfluorotetradecanoic acid	376-06-7		x
PFAS-PFMPA	Perfluoro-3-methoxypropanoic acid	377-73-1	x	
PFAS-8:2FTS	1H,1H, 2H, 2H-Perfluorodecane sulfonic acid	39108-34-4	x	
PFAS-PFTrDA	Perfluorotridecanoic Acid	72629-94-8		x
PFAS-9CI-PF1	9-chlorohexadecafluoro-3-oxanonane-1-sulfonic acid	756426-58-1	x	x
PFAS-4:2FTS	1H,1H, 2H, 2H-Perfluorohexane sulfonic acid	757124-72-4	x	
PFAS-11CI-PF2	11-chloroeicosafluoro-3-oxaundecane-1-sulfonic acid	763051-92-9	x	x
PFAS-PFMBA	Perfluoro-4-methoxybutanoic acid	863090-89-5	x	
PFAS-ADONA	4,8-dioxa-3H-perfluorononanoic acid	919005-14-4	x	x



SPEX.COM

+1.732.549.7144

+1.800.LAB.SPEX