



QC Samples for Water Pollution

Our WP Quality Control Samples are manufactured and tested for homogeneity following the requirements set for proficiency testing by NELAC, NIST/NVLAP and EPA. We recommend you run QC's throughout the year to quickly identify and correct any problems to ensure greater confidence in your final PT results.

All inorganic samples are supplied as 21 mL ampules unless otherwise specified. All **organic** samples are supplied as 2.2 mL ampules unless otherwise specified.

WP ORGANIC

WP TOTAL ORGANIC HALIDE (TOX)

Total Organic Halides 100-1000 µg/L

QC Part#: QC-TOX-WP Price: \$55.00

WP VOLATILES (1.5 mL)

Acetone	5-200 µg/L
Acetonitrile	5-200 µg/L
Acrolein	5-200 µg/L
Acrylonitrile	5-200 µg/L
Benzene	8-75 µg/L
Bromodichloromethane	7-75 µg/L
Bromoform	11-75 µg/L
Bromomethane	10-200 µg/L
2-Butanone (MEK)	5-200 µg/L
Carbon Disulfide	5-200 µg/L
Carbon Tetrachloride	9-75 µg/L
Chlorobenzene	11-75 µg/L
Chloroethane	10-200 µg/L
2-Chloroethylvinylether	5-200 µg/L
Chloroform	10-75 µg/L
Chloromethane	10-200 µg/L
Dibromochloromethane	9-75 µg/L
1,2-Dibromo-3-chloropropane (DBCP)	5-200 µg/L
1,2-Dibromoethane (EDB)	5-200 µg/L
Dibromomethane	5-200 µg/L
1,2 Dichlorobenzene	8-65 µg/L
1,3 Dichlorobenzene	8-60 µg/L
1,4 Dichlorobenzene	7-65 µg/L
Dichlorodifluoromethane	5-200 µg/L
1,1-Dichloroethane	5-200 µg/L
1,2 Dichloroethane	10-75 µg/L
1,1-Dichloroethene	10-200 µg/L
cis-1,2-Dichloroethene	5-200 µg/L
trans-1,2-Dichloroethene	10-200 µg/L
1,2-Dichloropropane	10-200 µg/L



WP VOLATILES cont.

cis-1,3-Dichloropropylene	5-200 µg/L
Ethylbenzene	8-75 µg/L
2-Hexanone	5-200 µg/L
Methylene Chloride	8-75 µg/L
4-Methyl-2-Pentanone (MIBK)	5-200 µg/L
Styrene	5-200 µg/L
1,1,1,2-Tetrachloroethane	5-200 µg/L
1,1,2,2-Tetrachloroethane	10-200 µg/L
Tetrachloroethene	9-75 µg/L
Toluene	7-65 µg/L
1,1,2-Trichloroethane	10-200 µg/L
1,2,3-Trichloropropane	5-200 µg/L
trans-1,3-Dichloropropene	25-250 µg/L
1,1,1-Trichloroethane	9-75 µg/L
Trichloroethene	7-75 µg/L
Trichlorofluoromethane	25-200 µg/L
Vinyl acetate	5-200 µg/L
Vinyl Chloride	10-200 µg/L
Xylenes, Total	10-500 µg/L

QC Part#: QC-VOA-WP Price: \$81.00

WP BASE/NEUTRALS

Acenaphthene	10-200 µg/L
Acenaphthylene	10-200 µg/L
Aniline	10-200 µg/L
Anthracene	30-200 µg/L
Benzidine	10-200 µg/L
Benzo(a)anthracene	30-200 µg/L
Benzo butyl phthalate	30-200 µg/L
Benzo(b)fluoranthene	30-200 µg/L
Benzo(g,h,i)perylene	10-200 µg/L
Benzo(k)fluoranthene	10-200 µg/L
Benzo(a)pyrene	10-200 µg/L
Benzyl alcohol	10-200 µg/L
4-Bromophenyl-phenylether	10-200 µg/L
bis(2-Chloroethoxy)methane	10-200 µg/L
bis(2-Chloroethyl)ether	10-200 µg/L
bis(2-Chloroisopropyl)ether	30-200 µg/L
bis(2-ethylhexyl)phthalate	30-200 µg/L
Carbazole	10-200 µg/L
4-Chloroaniline	10-200 µg/L
4-Chlorophenyl-phenylether	30-200 µg/L
1-Chloronaphthalene	10-200 µg/L
2-Chloronaphthalene	10-200 µg/L
Chrysene	30-200 µg/L



WP BASE/NEUTRALS cont.

Dibenzo(a,h)anthracene	10-200 µg/L
Dibenzofuran	10-200 µg/L
1,2-Dichlorobenzene	10-200 µg/L
1,3-Dichlorobenzene	10-200 µg/L
1,4-Dichlorobenzene	10-200 µg/L
3,3'-Dichlorobenzidine	10-200 µg/L
Di-n-octylphthalate	10-200 µg/L
Diethylphthalate	30-200 µg/L
Dimethylphthalate	30-200 µg/L
Di-n-butylphthalate	10-200 µg/L
2,4-Dinitrotoluene	10-200 µg/L
2,6-Dinitrotoluene	10-200 µg/L
Fluoranthene	10-200 µg/L
Fluorene	10-200 µg/L
Hexachlorobenzene	10-200 µg/L
Hexachlorobutadiene	10-200 µg/L
Hexachlorocyclopentadiene	10-200 µg/L
Hexachloroethane	30-200 µg/L
Indeno(1,2,3,cd)pyrene	30-200 µg/L
Isophorone	30-200 µg/L
1-Methylnaphthalene	10-200 µg/L
2-Methylnaphthalene	10-200 µg/L
Naphthalene	30-200 µg/L
2-Nitroaniline	10-200 µg/L
3-Nitroaniline	10-200 µg/L
4-Nitroaniline	10-200 µg/L
Nitrobenzene	10-200 µg/L
N-Nitrosodiethylamine	10-200 µg/L
N-Nitrosodimethylamine	10-200 µg/L
N-Nitroso-di-n-propylamine	10-200 µg/L
N-Nitrosodiphenylamine	10-200 µg/L
Phenanthrene	30-200 µg/L
Pyrene	10-200 µg/L
Pyridine	10-200 µg/L
1,2,4-Trichlorobenzene	10-200 µg/L

QC Part#: QC-BN-WP Price: \$84.00

WP ACIDS

Benzoic acid	10-200 µg/L
4-Chloro-3-methylphenol	10-200 µg/L
2-Chlorophenol	10-200 µg/L
2,4-Dichlorophenol	10-200 µg/L
2,6-Dichlorophenol	10-200 µg/L
2,4-Dimethylphenol	10-200 µg/L
2,4-Dinitrophenol	10-200 µg/L
2-Methylphenol	10-200 µg/L



WP ACIDS cont.

3-Methylphenol	10-200 µg/L
4-Methylphenol	10-200 µg/L
2-Methyl-4,6-Dinitrophenol	50-200 µg/L
2-Nitrophenol	10-200 µg/L
3-Nitrophenol	10-200 µg/L
4-Nitrophenol	10-200 µg/L
Phenol	10-200 µg/L
Pentachlorophenol	10-200 µg/L
2,4,5-Trichlorophenol	50-200 µg/L
2,4,6-Trichlorophenol	10-200 µg/L

QC Part#: QC-ACIDS-WP Price: \$62.00

WP PESTICIDES

Aldrin	0.1-3 µg/L
4,4'-DDD	0.1-7 µg/L
4,4'-DDE	0.1-4 µg/L
4,4'-DDT	0.1-8 µg/L
Dieldrin	0.1-5 µg/L
Heptachlor	0.15-4 µg/L
Heptachlor Epoxide (beta)	0.1-3 µg/L
Methoxychlor	0.5-25 µg/L
Alpha-BHC	2-15 µg/L
Beta-BHC	0.5-25 µg/L
Delta-BHC	7.5-75 µg/L
Gamma-BHC (Lindane)	3-25 µg/L
Endosulfan I	1-25 µg/L
Endosulfan II	20-75 µg/L
Endosulfan Sulfate	5-25 µg/L
Endrin	2-75 µg/L
Endrin aldehyde	0.5-15 µg/L
Endrin Ketone	0.5-15 µg/L

QC Part#: QC-PEST-WP Price: \$74.00

WP CHLORDANE

Chlordane (total)	1-13 µg/L
-------------------	-----------

QC Part#: QC-CHLR-WP Price: \$53.00

WP TOXAPHENE

Toxaphene (total)	2-10 µg/L
-------------------	-----------

QC Part#: QC-TXP-WP Price: \$53.00



WP PAHs

Acenaphthene	1-25 µg/L
Acenaphthylene	1-25 µg/L
Anthracene	1-25 µg/L
Benzo(a)anthracene	1-25 µg/L
Benzo(b)fluoranthene	1-25 µg/L
Benzo(k)fluoranthene	1-25 µg/L
Benzo(g,h,i)perylene	1-25 µg/L
Benzo(a)pyrene	1-25 µg/L
Chrysene	1-25 µg/L
Dibenz(a,h)anthracene	1-25 µg/L
Fluoranthene	1-25 µg/L
Fluorene	1-25 µg/L
Indeno(1,2,3-cd)pyrene	1-25 µg/L
Naphthalene	1-25 µg/L
Phenanthrene	1-25 µg/L
Pyrene	1-25 µg/L

QC Part#: QC-PAH-WP Price: \$62.00

WP PCBs IN WATER (2 AMPULES @ 2.2 mL)

Aroclor 1016/1242	3.8-13 µg/L
Aroclor 1221	1-15 µg/L
Aroclor 1232	1.4-4 µg/L
Aroclor 1248	1.5-5.5 µg/L
Aroclor 1254	1.7-5.5 µg/L
Aroclor 1260	1.6-5 µg/L

Each sample contains one specific aroclor

QC Part#: QC-PCBW-WP Price: \$57.00

WP PCBs IN OIL (2 AMPULES @ 1.5 g)

Aroclor 1016/1242	17-50 mg/kg
Aroclor 1254	16-50 mg/kg
Aroclor 1260	12-50 mg/kg

Each sample contains one specific aroclor

QC Part#: QC-PCBO-WP Price: \$90.00

WP HERBICIDES

2,4-D	1-10 µg/L
Dicamba	1-10 µg/L
2,4,5-T	1-10 µg/L
Silvex	1-10 µg/L

QC Part#: QC-HERB-WP Price: \$65.00